


EcoSac 5mm Tissue Retrieval System

Bolder Surgical is the exclusive U.S. distributor of the 60 mL and 80 mL EcoSac 5mm Tissue Retrieval System from Espiner Medical. This system easily fits down any 5mm cannula and provides durability, access, and ease of use.



Durability.

The EcoSac 5mm Tissue Retrieval System features a uniquely coated ripstop nylon that is resistant to tearing if punctured.

Access.

Easily fitting down any 5mm cannula, the EcoSac provides options to downsize ports and the 60ml and 80ml sac sizes will fully deploy during laparoscopic surgery in small spaces.


Ease of Use.

With the unique size and drawstring closure, the sac can be easily manipulated and retrieved with a standard 5mm grasper.

Clinical Benefits

Retrospective study of 11,475 patients showed appendix **retrieval bag use was associated with a 40% decrease in intra abdominal infection rates.**¹

Use of a retrieval bag can reduce risk of port site infection and intra-abdominal abscess, thereby **saving a hospital thousands of dollars** in complication related costs.²



“The use of appendix retrieval bags should be the standard of care during laparoscopic appendectomy.”¹

References

1. Adam C. Fields, Pamela Lu, Deanna L. Palenzuela, Ronald Bleday, Joel E. Goldberg, Jennifer Irani, Jennifer S. Davids, Nelya Melnitchouk. Does retrieval bag use during laparoscopic appendectomy reduce postoperative infection? *Surgery*, Vol. 165, Issue 5.
2. Lasheen A E, Abd Elaziz O, Elaal S A, Alkilany M, Sieda B, et al. Surgical Wound Infections After Laparoscopic Appendectomy With or Without Using Reusable Retrieval Bag: A Retrospective Study, *J Minim Invasive Surg Sci*. 2016 ; 5(2):e36894. doi: 10.17795/minsurgery-36894.

Ordering Information

Espiner Medical Eco-Sac 5mm 60ml	Catalogue No.	Package	Insertion Rod	Volume	Shaft Diameter
	EMP60ECO-R-10	Case of 10 Bags	45 cm	60 ml	5 mm
Espiner Medical Eco-Sac 5mm 80ml	Catalogue No.	Package	Insertion Rod	Volume	Shaft Diameter
	EMP80ECO-R-10	Case of 10 Bags	45 cm	80 ml	5 mm

At Bolder Surgical, we strive to elevate expectations in surgery. Leveraging our expertise in microsurgical devices and energy delivery, we innovate to give surgeons options they’ve never had, making them more proficient while minimizing the surgical impact on every patient.

Elevate your expectations with durable surgical instruments like the EcoSac 5mm Tissue Retrieval System that allows surgeons to utilize smaller ports and work efficiently in small spaces.



331 S. 104th Street Suite 200, Louisville, CO 80027
 P: 720-287-7133 F: 720-287-7135
www.boldersurgical.com